




User's Manual

Electronic Personal Dosimeter NRF50



Introduction

This User's Manual is a document to explain the operation of the Electronic Personal Dosimeter NRF50 sample unit briefly. It is organized to provide descriptions of parts, functions and operational instructions for optimal use. Please make sure that you read this manual carefully before operation. If there are some requirements or improvements about this sample unit, please contact Fuji Electric representative.

Also, in the event of unit malfunction, contact Fuji Electric representative immediately.

Handling Precaution

Please read the following handling precautions to ensure that you use the sample unit safely and avoid injury/damages. Please read this User's Manual carefully to understand all the precautions before using the sample unit.

	Measures of Precaution
 Attention	<ul style="list-style-type: none"> • The Dosimeter is a precision instrument; do not drop it or subject it to impact. • Keep the Dosimeter in a plastic bag for protection when use in an environment where chemical fumes, splashes/steam, full of dust and wastes are present. • Handle the Dosimeter with clean, dry hands. If becomes tainted, clean it with dry cloth. • Do not place the Dosimeter and metal objects in the same pocket. It may cause the Dosimeter breaking. • Avoid use where high frequency noise. Pay attention when use near the following devices: <ol style="list-style-type: none"> 1. Mobile phone 2. Premises/local wireless phone such as <PHS> 3. High power transceiver, or like kind 4. Microwave oven 5. Radar 6. Welding machine 7. Any other spark discharging or high intensity radio wave emitting devices <p>Especially keep the Dosimeter at least 5cm away from any mobile/wireless phones</p> • When the battery level is critically low, read the displayed value within 10 minutes. • Use AA alkaline battery only. During replacement, align the battery polarities correctly. • Prior to disposal of the used battery, protect exposed terminals with insulating tape to prevent shorting that may cause possible heating, rupture, or burning. Otherwise, injury or fire may result. • Do not throw the Dosimeter or battery into a fire. Do not disassemble them. • Keep distance between the buzzer and ears to avoid the injury. (Buzzer makes the sound over 90dB) • Do not use the Dosimeter as a survey meter. • If a hard impact operating on the Dosimeter, it is fear of having a crack to the Dosimeter. In this case, it is capable of deteriorations to the capabilities in a waterproof and a resistance of radio wave.

Contents

1.	Overview	5
2.	Contents.....	6
2.1	Sample unit package	6
2.2	Model.....	6
3.	Precaution.....	7
3.1	Operational conditions	7
3.2	Other requirements	7
4.	Description of Parts and Functions	8
4.1	Part names.....	8
4.2	Display function.....	10
4.3	Buzzer function	11
5.	Parts Replacement	15
5.1	Battery replacement.....	15
6.	Operational Instruction	16
6.1	When starting to use	16
6.2	During use (Normal operation)	18
6.3	During use (When alarm is generated)	26
6.4	During use (When communicating)	28
6.5	After use.....	28
7.	Care and Maintenance	29
7.1	Daily check and maintenance items	29
7.2	Consumable supplies.....	29
8.	Specification	30
8.1	General Specification.....	30
8.2	Storage data.....	31
9.	Appendix	32
9.1	Trouble shooting table	32
9.2	Disposal.....	33
9.3	Calibration	34
9.4	Block diagram	36

1. Overview

The Electronic Personal Dosimeter NRF50 (hereinafter referred to as NRF50) is designed to provide measurement of personal dose equivalent of external exposure to radiations (hereinafter referred to as dose).

NRF50 indicates accumulated dose or dose rate. If measured dose (rate) value exceeds preset dose (rate) threshold (Alarm threshold), NRF50 will be audibly alarmed and flashes LED.

Using the Setting Device and a PC, it is able to write PC-edited values to NRF 50 and read measurement trend data from the NRF50 via communication with the device.

If worn tight to the body, energy characteristic of the NRF50 enables direct reading of personal dose equivalent $H_p(10)$. (unit: rem)

2. Contents

2.1 Sample unit package

- | | |
|---------------------------------|---|
| (1) NRF50 | 1 |
| (2) Accessory | |
| ▪ Battery (AA alkaline battery) | 2 |

2.2 Model

NRF50

3. Precaution

Attention

3.1 Operational conditions


Item	Conditions
Temperature range	-10 °C to +50 °C
Relative humidity	95 % or less (No condensation)
Storage temperature	-25 °C to +50 °C

3.2 Other requirements

- (1) See User's Manual of "Setting Device" for information on parameter writing and data reading via the device and a PC.
- (2) Try turn OFF & ON the NRF50 if you encounter technical problems. See the "Troubleshooting Table" if not recovered.
- (3) Accumulated dose cannot be reset if the power-on-reset setting of NRF50 is "OFF". In this case, please reset the accumulated dose through the Setting Device.

4. Description of Parts and Functions


4.1 Part names



Functions

1. Buzzer : Sounds when starting, alarm and monitoring sound is generated.
2. Flash LED : Red LED on top of LCD flashes during alarm generation.
3. Auxiliary LED : Green LED flashes along to monitoring sound. Red LED flashes along to alarm. Yellow LED flashes along to preliminary alarm.
4. LCD : Liquid Crystal Display Indicator.
5. Battery Cap : Cap of battery compartment. Tighten with screw.
6. Operating Button : Switch indication screen. Also, alarm threshold can be changed by this button.
7. Calibration Mark : Shows location of calibration point. (Sensor position)
8. USB Connector : For power supply or communication by connecting USB cable.
9. IR : For communication with infrared communication device.
10. Call Button : Emergency alarm button. Press for more than 3 seconds to generate warning sound and indication.
11. EPD Information : For indicating model and serial number of the NRF50.


4.2 Display function



4.3 Buzzer function

4.3.1 Audible signal

Audible signals sound under the following circumstances:

Signal beep	Starts when	Beep pattern
Turn ON/OFF	NRF50 is turned on	
	NRF50 is turned off	No sound
Contact to setting device	Communication starts.	No sound
Data transmission	Successful completion of data transmission	No sound
	Setting values of NRF50 are changed using a setting device.	No sound
	Data transmission failed	No sound

4.3.2 Audible alarms

Alarm activations and beep patterns are as follows:

① Dose value alarm

Alarm type	LCD indication	Buzzer	Vibration	LED
Alarm for accumulated dose		Sounds 3 times per 1 second	Vibrates 1 time per 1 second	Flashes 1 time per 1 second
Preliminary alarm for accumulated dose				Alarm: Flash LED + Auxiliary LED (RED)
Alarm for dose rate		Sounds 2 times per 1 second (Long beep)		Preliminary alarm: Auxiliary LED (YELLOW)
Preliminary alarm for dose rate				
<p>Operating time</p>				

② Operation time alarm

Alarm type	LCD indication	Buzzer	Vibration	LED
Time alarm		Sounds 1 time per 1 second (Short sound)	Vibrates 1 time per 10 seconds	Flashes 1 time per 10 seconds
<p>Operating time</p>				

③ Low battery voltage

Alarm type	LCD indication	Buzzer	Vibration	LED
Low battery voltage		Sounds 3 times per 10 minutes (Short sound)	No vibration	Flashes 1 time per 10 minutes Auxiliary LED (RED)
Operating time 				

④ Indication of abnormality

Alarm type	LCD indication	Buzzer	Vibration	LED
Detector optical error		Sounds 4 times per 1 second	Vibrates 1 time per 1 second	Flashes 1 time per 1 second
Memory error				Flash LED + Auxiliary LED (RED)
Memory error (When measurement stops)				
RTC error				
RTC error (When measurement stops)				
Emergency alarm				
Operating time				

4.3.3 Monitoring Sound

Buzzer sounds 1 time, if accumulated dose reaches a preset value of dose interval for monitoring sound. Monitoring sound interval can be chosen from 6 types such as “OFF”, “10mrem”, “1mrem”, “0.1mrem”, “8 counts” and “4 counts”. Please see User’s Manual of “Setting Device” for details.

5. Parts Replacement

5.1 Battery replacement

Follow these steps to replace the batteries:

(1) Press and hold “Ⓢ” button to turn NRF50 off.

(2) Open the cap of battery compartment using screwdriver.


(3) Replace the batteries.

Insert new ones properly in the compartment with care to check the directions “+” and “-”.

(Both batteries are positioned in the same direction)

(4) Close the cap and tighten with screwdriver.






 Attention	<ul style="list-style-type: none">① When replacing batteries, make sure to turn off NRF50.② During replacement, align the battery polarity correctly.③ Use only AA Alkaline battery.
---	--

6. Operational Instruction

6.1 When starting to use

- (1) Press and hold “⊙” button for more than 3 seconds to start the NRF50. Confirm the power is ON (one beep) and LCD.



Check items	Confirmation method	
Audible signal (1 beep)	Generates when turning on the power.	
Indicated dose value	0.0 mrem or 0.0 mrem/h (Accumulated dose may not be 0.0 mrem if the power-on-reset setting of NRF50 is “OFF”.)	
LCD	<p>[Normal display]</p>  <p>When accumulated dose is 0.0mrem</p>	<p>[Abnormal display]</p>  <p>When generating alarm for low battery</p>


(2) White backlight turns on by selecting any of “◎”, “△” or “▽” button.

Display of accumulated dose (rem) and dose rate (rem/h) can be switched by pushing the button of “△” or “▽”. “◎” button can also switch the display to other screens.


Please see section 6.2 in details for screen change methods by “◎” button.

(3) Put NRF50 in the chest pocket as shown below.

① When putting in an outside pocket



② When putting in an








※Direction of NRF50:




As EPD is viewed from wearer, operation button and auxiliary LED must be positioned so that they are in the right side, and buzzer faces outward.



6.2 During use (Normal operation)

LCD display during operation of button

※LCD display will be changed when pressing the designated button with backlight flashed.

No.	Display	Remarks
1	<p>【Basic display】 ⇒To No.2</p>  <p>【Number indicating display】 ⇒To No.3</p> <p>Short press of “◎”</p>  <p>【Setting display】 ⇒To No.4</p> <p>Short press of “◎”</p>  <p>Short pressing of “◎”</p>	<p>Power turns off by long press of “◎” button. (Basic display or number indicating display)</p> <p>Indication is changed by short press of “△” or “▽” button in each display.</p>
2	<p>【Basic display】 Switching</p>  <p>Indication of accumulated dose</p> <p>Short press of “△” or “▽”</p>  <p>Indication of dose rate</p>	<p>Power turns off by long press of “◎” button.</p>





No.	Display	Remarks
3	<p data-bbox="240 286 745 322">【Number indicating display】 Switching</p> <p data-bbox="228 356 459 392">Short press of “▽”</p> <div data-bbox="311 423 681 586">  <p data-bbox="707 562 997 593">EPD number indication</p> </div> <p data-bbox="228 627 459 663">Short press of “▽”</p> <p data-bbox="544 627 775 663">Short press of “△”</p> <div data-bbox="311 694 681 857">  <p data-bbox="707 833 967 864">ID number indication</p> </div> <p data-bbox="228 898 459 934">Short press of “▽”</p> <p data-bbox="544 898 791 934">Short press of “△”</p> <div data-bbox="311 965 681 1128">  <p data-bbox="707 1104 1003 1135">RWP number indication</p> </div> <p data-bbox="552 1200 783 1236">Short press of “△”</p>	<p data-bbox="1193 286 1485 353">Power turns off by long press of “◎” button.</p>

No.	Display	Remarks
4	<p data-bbox="240 286 596 322">【Setting Display】 Switching</p> <p data-bbox="228 356 459 392">Short press of “▽”</p> <div data-bbox="312 423 683 586">  </div> <p data-bbox="708 443 1032 573">Configuration Long press of “◎” to change display for setting items</p> <p data-bbox="228 694 459 730">Short press of “▽”</p> <p data-bbox="571 694 971 730">Short press of “△” ⇒ To No.5</p> <div data-bbox="312 792 683 956">  </div> <p data-bbox="708 815 1040 945">Memo Long press of “◎” to change a memo indicating display</p> <p data-bbox="564 999 968 1034">Short press of “△” ⇒ To No.6</p>	<p data-bbox="1193 286 1458 356">Long press of “◎” to change each display.</p>

No.	Display	Remarks
5	<p data-bbox="240 286 667 320">【Setting items display】 Switching</p> <p data-bbox="240 353 491 421">Short press of “▽”</p> <div data-bbox="311 434 679 595"> <p>d 0.0 Hp10 NAME mrem JOHN MAYER</p> </div> <p data-bbox="719 499 884 533">Name setting</p> <p data-bbox="240 640 491 707">Short press of “▽”</p> <p data-bbox="560 640 788 674">Short press of “△”</p> <div data-bbox="311 719 679 880"> <p>d 0.0 Hp10 DOSE mrem 00. 1000 rem</p> </div> <p data-bbox="719 770 948 837">Accumulated dose alarm setting</p> <p data-bbox="240 922 491 990">Short press of “▽”</p> <p data-bbox="560 922 788 956">Short press of “△”</p> <div data-bbox="311 1003 679 1164"> <p>d 0.0 Hp10 DOSE mrem 00. 0500 rem</p> </div> <p data-bbox="719 1048 1027 1115">Accumulated dose preliminary alarm setting</p> <p data-bbox="240 1205 491 1272">Short press of “▽”</p> <p data-bbox="560 1205 788 1238">Short press of “△”</p> <div data-bbox="311 1285 679 1447"> <p>d 0.0 Hp10 RATE mrem 00. 0100 rem/h</p> </div> <p data-bbox="719 1330 884 1397">Dose rate alarm setting</p> <p data-bbox="240 1489 491 1556">Short press of “▽”</p> <p data-bbox="560 1489 788 1523">Short press of “△”</p> <div data-bbox="311 1570 679 1731"> <p>d 0.0 Hp10 RATE mrem 00. 0050 rem/h</p> </div> <p data-bbox="719 1637 1027 1704">Dose rate preliminary alarm setting</p> <p data-bbox="560 1778 788 1812">Short press of “△”</p> <p data-bbox="427 1845 628 1879">to the next page</p>	<p data-bbox="1187 286 1474 421">Short press of “◎” to change set values. (Yellow backlight turns on)</p> <p data-bbox="1187 454 1442 521">Long press of “◎” to confirm the change.</p> <p data-bbox="1187 555 1458 645">Name can be entered with 10 characters from A to Z.</p> <p data-bbox="1187 678 1474 745">Alarm value can be set with numbers of 0 to 9.</p>

No.	Display	Remarks
5	<p>【Setting items display】 Switching (continued)</p> <p>from the previous page</p> <p>Short press of “▽”</p> <div data-bbox="312 499 681 660"> <p>d 0.0 Hp10 TIME ALARM mrem 08:00</p> </div> <p>Operating time alarm setting</p> <p>Short press of “▽” Short press of “△”</p> <div data-bbox="312 784 681 945"> <p>d 0.0 Hp10 BEEP LEVEL mrem 10mrem</p> </div> <p>Monitoring sound setting</p> <p>Short press of “▽” Short press of “△”</p> <div data-bbox="312 1068 687 1232"> <p>d 0.0 Hp10 DATE SETUP mrem 08/ 20/ 2014</p> </div> <p>Date setting</p> <p>Short press of “▽” Short press of “△”</p> <div data-bbox="312 1355 687 1520"> <p>d 0.0 Hp10 TIME SETUP mrem 10: 35: 00</p> </div> <p>Time setting</p> <p>Short press of “▽” Short press of “△”</p> <div data-bbox="312 1644 681 1807"> <p>d 0.0 Hp10 MODE STANDALONE mrem</p> </div> <p>Mode setting (Only confirmation)</p> <p>Short press of “△”</p> <p>to the next page</p>	<p>Change a set value for each items by short press of “◎”. (Yellow backlight turns on)</p> <p>Confirm the change by long press of “◎”.</p> <p>Alarm setting, date and time setting can be set with numbers 0 to 9.</p> <p>Monitoring sound can be set in “10mrem”, “1mrem”, “0.1mrem”, “8PLS” (8 counts), “4PLS” (4 counts) and “OFF”.</p>







No.	Display	Remarks
5	<p>【Setting items display】 Switching (continued)</p> <p>from the previous page</p> <p>Short press of “▽” ↓ ↑ Short press of “△”</p> <div data-bbox="312 501 681 660"> </div> <p>Version (Only confirmation)</p> <p>Short press of “▽” ↓ ↑ Short press of “△”</p> <div data-bbox="312 784 681 943"> </div> <p>Setting completed ⇒ To No.4</p> <p>Short press of “△” ↑</p>	<p>Return to setting display by short press of “◎”.</p>







No.	Display	Remarks
6	<p data-bbox="236 286 443 320">【Memorandum】</p> <div data-bbox="308 353 675 517">  </div> <p data-bbox="707 409 954 477">Memo display (1 to 10 characters)</p> <div data-bbox="236 528 786 618"> <p>Short press of “▽” ↓ ↑ Short press of “△”</p> </div> <div data-bbox="308 629 675 792">  </div> <p data-bbox="707 685 970 752">Memo display (11 to 20 characters)</p> <div data-bbox="236 808 786 898"> <p>Short press of “▽” ↓ ↑ Short press of “△”</p> </div> <div data-bbox="308 909 675 1072">  </div> <p data-bbox="707 954 970 1021">Memo display (21 to 30 characters)</p> <div data-bbox="236 1088 786 1178"> <p>Short press of “▽” ↓ ↑ Short press of “△”</p> </div> <div data-bbox="308 1189 675 1352">  </div> <p data-bbox="707 1234 970 1301">Memo display (31 to 32 characters)</p>	<p data-bbox="1185 286 1509 353">Return to setting display by long press of “◎”.</p>

6.3 During use (When alarm is generated)

LCD display when generating alarm





※See Chapter 4 for operation of buzzer, vibration and LED during alarm generation.

No	Item	Display	Remarks
1	Accumulated dose alarm	<p>Alarm</p>  <p>Preliminary Alarm</p> 	Preliminary alarm is generated, when exceeding warning set value. Alarm is generated, when exceeding alarm set value.
2	Dose rate alarm	<p>Alarm</p>  <p>Preliminary alarm</p> 	Preliminary alarm is generated, when exceeding warning set value. Alarm is generated, when exceeding alarm set value. See chapter 4 for alarm operation. Alarm is canceled, when dose rate reaches 90% of preset value or less.
3	Operation time alarm		When exceeding alarm set time, time alarm is generated.
4	Low battery voltage		Low battery voltage alarm is generated, when battery voltage reaches less than 1.1V. Remaining operational hours is displayed with numbers 1 to 9 h, then power OFF after 1 hour.

No	Item	Display	Remarks
5	Detector optical error	 <p>The display shows a red background with a radiation symbol in the top left. At the top, it displays the date '08/20/2014' and time '10:24'. Below this, it shows 'R' followed by five dashes '-----'. In the top right corner, there is a battery level indicator and the text 'Hp10 rem/h'. The main display area shows 'OPTI ERROR' in large, bold, black letters.</p>	When detector cannot work properly for internal optical pulse check by internal LED, detector optical error alarm is generated.
6	Memory error	<p>When not affecting counting</p>  <p>The display shows a red background with a radiation symbol in the top left. At the top, it displays the date '08/20/2014' and time '10:24'. Below this, it shows 'R' followed by the value '0.300'. In the top right corner, there is a battery level indicator and the text 'Hp10 rem/h'. The main display area shows 'MEM ERROR' in large, bold, black letters.</p> <p>When affecting counting</p>  <p>The display shows a red background with a radiation symbol in the top left. At the top, it displays the date '08/20/2014' and time '10:24'. Below this, it shows 'R' followed by five dashes '-----'. In the top right corner, there is a battery level indicator and the text 'Hp10 rem/h'. The main display area shows 'MEM ERROR' in large, bold, black letters.</p>	When memory failure occurs during data backup, memory error alarm is generated. When not affecting counting, operation continues, and when affecting counting, operation stops holding the dose (rate) display.
7	RTC error	<p>When not affecting counting</p>  <p>The display shows a red background with a radiation symbol in the top left. At the top, it displays the date '08/20/2014' and time '10:24'. Below this, it shows 'R' followed by the value '0.101'. In the top right corner, there is a battery level indicator and the text 'Hp10 rem/h'. The main display area shows 'RTC ERROR' in large, bold, black letters.</p> <p>When affecting counting</p>  <p>The display shows a red background with a radiation symbol in the top left. At the top, it displays the date '08/20/2014' and time '10:24'. Below this, it shows 'R' followed by five dashes '-----'. In the top right corner, there is a battery level indicator and the text 'Hp10 rem/h'. The main display area shows 'RTC ERROR' in large, bold, black letters.</p>	RTC error alarm is generated, Real Time Clock (RTC) IC error occurs. When not affecting counting, operation continues, and when affecting counting, operation stops holding the dose (rate) display.
8	Emergency alarm	 <p>The display shows a red background with a radiation symbol in the top left. At the top, it displays the date '08/20/2014' and time '10:24'. Below this, it shows 'R' followed by the value '12.1'. In the top right corner, there is a battery level indicator and the text 'Hp10 mrem/h'. The main display area shows 'EMERGENCY' in large, bold, black letters.</p>	Emergency alarm is generated, when pressing "call" button for more than 3 seconds.

6.4 During use (When communicating)

LCD display during communication

No	Item	Display	Remarks
1	Infrared communication		"IR" is displayed on the upper left of screen during infrared communication.
2	USB communication		Icon does not appear during USB communication.
3	Bluetooth communication		"Bluetooth mark" is displayed on the upper left of screen during Bluetooth communication.
4	900MHz wireless telemetry		"Telemetry mark" is displayed on the upper left of screen during communication with 900MHz wireless telemetry.

6.5 After use

Turn OFF by long press of "Ⓢ" on basic display or number indicating display.

7. Care and Maintenance

Check the NRF50 as specified below to ensure quality performance.

7.1 Daily check and maintenance items

No	Check items	Procedures	Check point
1	Appearance	Check the NRF50 visually. When to check ; Before use and after battery replacement Check purpose ; Check if there is no abnormality with the case and battery cap.	No signs of crack, damages or breakage on the case. No signs of gap between case and battery cap.
2	Indication error/calibration	To confirm the indication error within 10% to the reference dose equivalent using Cs-137. When to check: 1 year or less Check purpose: To optimize the dose management.	If indication error is over 10%, please Contact Fuji Electric representative for calibration of NRF50.

7.2 Consumable supplies

Please Contact Fuji Electric representative for the following consumables.

- ① Battery cap with O-ring : Order from 1 unit
- ② USB connector cap : Order from 1 unit

8. Specification

8.1 General Specification

Model : NRF50
Category : Grade 1
Detector : Silicon semi-conductor
Radiation type : γ (X) rays (30 keV to 7.0 MeV)
Dose display range: 0.1 mrem to 1000 rem, 0.1 mrem/h to 1000 rem/h
Effective measurement range: 2.0 mrem to 1000 rem, 0.05 mrem/h to 1000 rem/h (accumulated dose)
10.0 mrem/h to 1000 rem/h (dose rate)
Rated range : Gmha (category for IEC61526)
60 keV to 6 MeV, 10.0 mrem to 1000 rem, 0.05 mrem/h to 1000 rem/h
Indication error : Within $\pm 10\%$ (Cs-137, 2.0 mrem to 1000 rem)
Within $\pm 30\%$ (Cs-137, 10.0 mrem/h to 100.0 mrem/h)
Within $\pm 20\%$ (Cs-137, 100.0 mrem/h to 1000 rem/h)
Energy response : Relative response 0.71 to 1.67 (60 keV to 6 MeV, 0° to 60°)
Temperature characteristic: -13 % to +18 % (20 °C, -10°C to +40°C)

Temperature test result

Left for 4 hours in each temperatures. Evaluated by amount of change in the last 30 minutes.

Reference (20°C)		-10°C			50°C		
variation	(mrem)	variation	(mrem)	result	variation	(mrem)	result
	1.004		1.032	2.79%		0.932	-7.17%

Dimensions : approx. 105 mm(H) × 60 mm(W) × 29/18 mm(D) (excluding protrusion)
Weight : approx. 170 g (2 Batteries), approx. 100 g (Battery excluded)
Battery : AA alkaline battery (x 2)
Continuous operating time: more than 1000 hours (under normal temperature, no alarms, new battery)
Reference standards: IEC61526 Ed3.0(2010), ANSI N42.20(2003)

8.2 Storage data

① List of storage data (Updated value is stored in EEPROM every 1 minute)

- EPD number
- Current time
- Current accumulated dose
- Current dose rate
- Operating time
- Alarm setting values (Accumulated dose, Dose rate : 2 for each)
- Time alarm setting value
- Calibration factor
- Error flag
- Condition flag
- Other setting values

② Trend data storage

Following data is stored at preset interval of trend. (max. 4000)

- Measurement date
- Accumulated dose
- Maximum dose rate

9. Appendix

9.1 Trouble shooting table

Error Indication	Possible Cause	Suggested Solution
「OPTI ERROR」	(1)Sensor unit malfunction (2)CPU malfunction	(1),(2) Contact Fuji Electric representative.
「MEM ERROR」	(1)EEPROM malfunction (2)CPU malfunction	(1),(2) Contact Fuji Electric representative.
「RTC ERROR」	(1)RTC malfunction (2)CPU malfunction	(1),(2) Contact Fuji Electric representative.

When returning the item to Fuji Electric representative, please provide with precise details of problems.

Note : This table is applied only to the malfunctions that occurs during use.

Symptom	Possible Cause	Suggested Solution
No indications on LCD	(1) Defective battery connection (2) Mode switch malfunction (3) LCD malfunction (4) CPU malfunction	(1) Check battery polarity and there is no exogenous material. (2)~(4) Contact Fuji Electric representative.
Characters on LCD are garbled.	(1) LCD malfunction (2) CPU malfunction	(1),(2) Contact Fuji Electric representative.
Backlight does not light when pressing a mode switch.	(1) Mode switch malfunction (2) LCD malfunction (3) CPU malfunction	(1)~(3) Contact Fuji Electric representative.
Dose error • Dose accumulation does not work • Displayed dose is high • Displayed dose is low	(1) LCD malfunction (2) Sensor unit malfunction (3) CPU malfunction (4) Calibration constant trouble	(1)~(3) Contact Fuji Electric representative. (4) Check calibration factor. Contact Fuji Electric representative for calibration method.
Buzzer does not sound	(If display is correct) (1) Attached Exogenous material (2) Set frequency failure (3) Buzzer lead wire is broken (4) CPU malfunction	(1)Check if there is no dust. If it is not improved after removing a dust, Contact Fuji Electric representative. (2)~(4) Contact Fuji Electric representative.
Vibration does not work	(1) Vibration malfunction (2) CPU malfunction	(1),(2) Contact Fuji Electric representative.

Symptom	Possible Cause	Suggested Solution
LED does not light	(1) LED unit malfunction (2) CPU malfunction	(1),(2) Contact Fuji Electric representative.
Operational hour is short Battery voltage alarm is always displayed	(1) End of the battery's life (2) Increase the current consumption (3) CPU malfunction (4) Failure of voltage decline detection value	(1) Replace with a new battery. See "5.1". (2) Check the proper contacts in the battery compartment and there is no exogenous material in the battery case. If trouble continues, Contact Fuji Electric representative. (3),(4) Contact Fuji Electric representative.
IR communication is unable.	(1) Communication distance is too far (2) Communication port is dirty. (3) CPU malfunction (4) Malfunction of setting software (PC)	(1) Set the distance between Communication port of NRF50 and the Setting device within 5cm. Also confirm that these windows are face to face. (2) Clean the communication part with soft cloth. Check if there is no exogenous material. If trouble continues, contact Fuji Electric representative. (3) Contact Fuji Electric representative. (4) Refer to the instruction manual of setting software to check if setting is proper.
USB communication does not work	(1) Cable malfunction (2) Connector malfunction (3) CPU malfunction	(1) Check conduction to confirm there is no problem with the cable. (2),(3) Contact Fuji Electric representative.
Bluetooth communication does not work	(1) Module malfunction (2) CPU malfunction (3) Abnormal pairing setting	(1),(2) Contact Fuji Electric representative. (3) Refer to software instruction manual to check if setting is proper.
Telemetry communication does not work	(1) Module malfunction (2) CPU malfunction	(1),(2) Contact Fuji Electric representative.
Crack, breakage, damage on the case	(1) Breakage due to drops, etc.	(1) Contact Fuji Electric representative.

9.2 Disposal

Please follow the local law and regulation for disposal of the product.

9.3 Calibration

This section describes the calibration procedures for NRF50.

Expose NRF50 to the gamma-ray sources such as ^{137}Cs and ^{60}Co .


A dose should be measured by placing the source at a certain distance (calibration distance) from reference point of NRF50 so that true value of the dose is traceable to the National Standard.

(1) Determination of a reference dose rate value (R_0)

- Determine a reference dose rate value (R_0) by the following method:
 - a. Calculate R_0 from the reference source activity and the distance between the reference source and reference point of NRF50 (calibration distance).
 - b. Or the dose value at the reference point (R_0) may be simply well-known by field calibration/characterization.

(2) Dose value (R_1) measurement

- Place the source such as ^{137}Cs and ^{60}Co at the calibration distance from reference point of NRF50.
- Take the dose reading (R_1) after irradiation which gives enough value of dose.



Example of Geometrical Conditions

(3) Calculation of the calibration factor

- Compare the reference dose (R_0) and the dose reading (R_1). If there is an unacceptable difference between R_0 and R_1 , change the calibration factor.

In general, the calibration factor (C_1) is calculated by the following formula:


$$C_1 = C_0 \times R_0/R_1$$

C_0 : Original Calibration Factor

(4) Setup of the calibration factor

- To change the calibration factor, perform the following procedures:
 - a. After the irradiation, run configuration software.
 - b. Click on "Calibration", enter the calculated calibration factor (C_1) to the new value of gamma-ray calibration factor.
 - c. Press the "Write" button.
 - d. Confirm the current value is set to the new value.

9.4 Block diagram





★ Your Comment ★

Dear customers,

Any comments/ requests/ suggestions regarding our instruction manual? Please feel free to contact us by filling out this form and give to our sales representative.

Document No	TN5A2293	Date	year month day		
Name of manual	Electronic Personal Dosimeter NRF50 User's Manual	Submitted by	Company		
			Dept		
			Name		

Page	Line	Description
		Comments, requests, questions, doubts (error, more info needed, terminology not consistent, typos)Put with O

Manufacturer	Received by		Received date		year month day	Receipt #	
--------------	-------------	--	---------------	--	----------------	-----------	--